


# Crypto-Strip

Rapid diagnostic test for *in vitro* detection of *Cryptosporidium* in stool specimens

- Easy to use
- Easy to read
- Specific
- Sensitive
- Reliable
- Rapid format (10')



## Accurate, economical and fast procedure!



## Performance

compared to EIA (n=103)

compared to immuno-enzymatic rapid test (n=100)

<b>Specificity</b>	<b>100% (85/85)</b>
<b>Sensitivity</b>	<b>94.4% (17/18)</b>
<b>PPV</b>	<b>100% (17/17)</b>
<b>NPV</b>	<b>98.8% (85/86)</b>
<b>Accuracy</b>	<b>99% (102/103)</b>

<b>100% (53/53)</b>
<b>95.7% (45/47)</b>
<b>100% (45/45)</b>
<b>96.4% (53/55)</b>
<b>98% (98/100)</b>

## Bibliography

For more information:

[www.corisbio.com/public/science](http://www.corisbio.com/public/science)

Or send an e-mail to: [info@corisbio.com](mailto:info@corisbio.com)

## Ordering Information

### Enteric range (25 tests)

Description	Code
Rota-Strip	C-1001
Adeno-Strip	C-1002
40/41 Adeno-Strip	C-1003
Combi-Strip	C-1004
Crypto-Strip	C-1005
O157 Coli-Strip (20 tests)	C-1011
Giardia-Strip	C-1013
GastroVir-Strip	C-1016
Crypto/Giardia Duo-Strip	C-1018

### Respiratory range (25 tests)

Description	Code
RSV Respi-Strip	C-1006
Adeno Respi-Strip	C-1009
Influ-A Respi-Strip	C-1010
Influ-A&B Respi-Strip	C-1012
Legionella V- Test*	C-1815

All our products are also available in packaging of 10 tests, and in packaging of 20 tests, individually packed.

\*Only available in 10 or 20 tests.



Science Park CREALYS  
Rue Jean Sonet 4A - 5032 Gembloux - Belgium  
Ph.: +32 (0)81 719 917 - Fx.: +32 (0)81 719 919  
[info@corisbio.com](mailto:info@corisbio.com) - [www.corisbio.com](http://www.corisbio.com)

Distributed by:

